# New top-lid packaging solution for rigid blister applications

# Polybond



#### CONTEXT

Over the years, multiple materials have been developed to serve as a potential cost effective alternative to Tyvek® based packaging. The first market developments have been observed in the field of flexible formable packaging where coated and surface treated webs made of a mix of cellulose and synthetic binders or fibers have progressively gained markets shared.

Among manufacturing devices that have been transferred, one can mention drapes, gowns, catheters & IV Sets...

At the end of the year 2000s, a first attempt had been made with the launch of Ovantex® by Oliver, which was a product development associating

- $\boldsymbol{\cdot}$  the coating technology of the by-then Oliver company
- to a versatile cellulose-based web material encompassing synthetic fibers & binders to enhance material mechanical strength & cohesion
- the total gauge of the material was of 101  ${
  m g/m^2}$
- enhanced mechanical properties
- material had been tested and complied with requirements for EO & Irradiation sterilization

In the early 2010s, polymer-reinforced webs were introduced with surface treatment technology helping to drive costs down by removing the need for additional coated layer on base web material.

The solution could be used on regular grades of film, granted that they were comprising a standard layer of Polyethylene for sealing & peeling performances, and has quickly found its market in the field of industrial peel pouch & roll stock packaging solutions for drapes, gowns, lap sponges and IV sets.

Nevertheless, specific requirements of sealing to rigid packaging could not be answered by this new solution, due to

- seal behavior to PVC, APET, PETG or HIPS being different from the one of a Pe- EVA
- peel angle and mechanical stress when opening a rigid blister being more challenging from aseptic opening requirements than the flexible packaging solutions.

#### NEW PRODUCT DEVELOPMENT

Hence the development over the last 3 years of new materials aiming at bridging the gap and delivering the right level of performance for blister needs:

- stronger intense cohesion of material
- flexibility and stretch enhancement for large peel angle requirements
- bespoke coated adhesives adapted to the different film natures that can be found on market

# POLYBOND AGP 90 gsm & POLYBOND BHP 85 gsm are the result of these product developments to satisfy unanswered customer needs.

**POLYBOND AGP 90 gsm** for sealing to PVC, which is still commonly used in Asia and on some European market, and to APET

2 POLYBOND BHP 85 gsm for sealing to HIPS



MATERIAL: fiber mesh of cellulose and polymer binders, with grid coating pattern



TECHNICAL PUBLICATION

## **Polybond**<sup>-</sup> AGP 90 gsm





### FOCUSING ON POLYBOND AGP 90 gsm TESTING & RESULTS



#### • Typical properties

**POLYBOND AGP 90 gsm** has been designed for sealing to PVC and APET rigid material and tested on several web thicknesses, from 300 to 700 μ, and form fill seal packaging lines. It's unique fiber mesh structure mixes cellulose together with a high content level of synthetic binders.

It allows for a very flexible material, and for superior performance in regard with prevention of the risk of break, puncture & burst. The choice of a water-based coating technology in addition to a significant share of renewable resource based raw material also make it a quite environmental friendly solution to your packaging needs. Available in roll stock or for converting into die cutted top lid.

#### • Illustrating material performance through packaging line testing

**POLYBOND AGP 90 gsm** is a coated reinforced web which complies with all requirements of international standards as illustrated by the below table detailing typical values and relevant test methods for the developed packaging material.

#### MATERIAL DESCRIPTION

|                     |                         |         | POROUS WEB                         |                                    |                                     |  |  |
|---------------------|-------------------------|---------|------------------------------------|------------------------------------|-------------------------------------|--|--|
| Product designation |                         |         | POLYBOND AGP 90 gsm                |                                    |                                     |  |  |
| Product description | I                       |         | GRID COATED POLYMER-REINFORCED WEB |                                    |                                     |  |  |
|                     | Standard Test<br>Method | Unit    | Typical value before sterilization | Typical value<br>EO* sterilization | Typical value<br>IR** sterilization |  |  |
| SUBSTANCE           | ISO 536                 | g/m²    | 90                                 | 90                                 | 90                                  |  |  |
| TENSILE STRENGTH    | ISO 1924-2              | kN/m    | 75                                 | 75                                 | 60                                  |  |  |
| BURST STRENGTH      | ISO 2758                | kPa     | 415                                | 385                                | 305                                 |  |  |
| TEAR STRENGTH       | ISO 1974                | mN      | 630                                | 650                                | 535                                 |  |  |
| COBB TEST           | ISO 535                 | g/m²    | 18.7                               | 17.0                               | 16.5                                |  |  |
| MAX PORE SIZE       | EN 868-6 app C          | μm      | 18.7                               | 19.7                               | 18.6                                |  |  |
| AIR PERMEANCE       | ISO 5636-3              | µm/Pa.s | 1.7                                | 2.3                                | 1.7                                 |  |  |
|                     |                         |         |                                    |                                    |                                     |  |  |

\* EO: Time 6h, Temperature 35°C at 50 % RH, EO concentration 1 Kg/m<sup>3</sup>. \*\* Irradiation : 28.6 kGy



# **Polybond**<sup>-</sup> AGP 90 gsm



Polybond AGP 90 gsm has been recently evaluated on:

- Multivac R 145 model
- To APET and PVC film materials
- In thicknesses of 600 and 300  $\mu$  respectively

60 different blisters have been manufactured and tested, opened & characterized. Mechanical, seal strength and 180° peel angle peel opening properties have been evaluated before and after sterilization. Sterilization conditions that have been used are as follow:

- EO on 6 hours with a temperature 35°C at 50 % RH and a concentration of 1 Kg/m<sup>3</sup>
- Irradiation with 28.6 kGy

Running conditions on packaging lines have been:

- Range of sealing temperature: 110-160°C 230-320 Fahrenheit
- Tested dwell time: 1,6-2,4 second
- Range of pressure: 300-500 KPa

Results of design of experiment:

- Medium to high seal strength levels have been achieved both when sealing to APET and PVC material
- Pressure applied on packaging line during manufacturing process has no noticeable impact on packaging seal and peel performance (real gas test diameter: 60 shore)
- A wide operating window (temperature, pressure, dwell time) has been identified as illustrated through the below charts

PVC 300µ Sealing strength (cN/15mm):



#### **DESIGNED FOR A LARGE OPERATING WINDOW ON YOUR PACKAGING LINES**

Design of Experiment of Polybond AGP 90 gsm to PVC 300µ sealing strength, 60 tested samples



Design of Experiment of Polybond AGP 90 gsm to APET  $600\mu$  sealing strength, 60 tested samples

Pressure (bars) / Temperature (°C) 5.0 SEALING STRENGTH 4.5 < 300 (bars) 300 - 360 > 360 Pressure 4,0 FIEXED Time 2s 3,5 3,0 120 130 140 110 Temperature (°C)

Design of Experiment of Polybond AGP 90 gsm to PVC 300  $\mu$  sealing strength, 60 tested samples



Design of Experiment of Polybond AGP 90 gsm to APET  $600\mu$  sealing strength, 60 tested samples



# **Polybond**<sup>-</sup> AGP 90 gsm



More data on product testing and results:

. . . . . . . . . . . . . . . . . . .

#### STERILE BARRIER SYSTEM PERFORMANCES

|                            |                                  |         | POROUS WEB                          |                                       |                                      |  |  |
|----------------------------|----------------------------------|---------|-------------------------------------|---------------------------------------|--------------------------------------|--|--|
| Product designation        |                                  |         | POLYBOND AGP 90 gsm                 |                                       |                                      |  |  |
| Product description        |                                  |         | GRID COATED POLYMER-REINFORCED WEB  |                                       |                                      |  |  |
|                            |                                  | Unit    | Before<br>sterilization<br>PVC 300µ | After EO<br>sterilization<br>PVC 300µ | Before<br>sterilization<br>APET 600µ | After Irrad.<br>sterilization<br>APET 600µ |  |
| TESTED CONDITIONS          | Tray Dimensions                  | mm      | 337x151x20                          | 337x151x20                            | 156x136x20                           | 156x136x20                                 |  |
| TYPICAL USES<br>CONDITIONS | Sealing temperature              | °C      | 130                                 | 130                                   | 140                                  | 140  |  |
|                            | Dwell time                       | S       | 2                                   | 2                                     | 2                                    | 2  |  |
|                            | Sealing pressure                 | KPa     | 400                                 | 400                                   | 400                                  | 400  |  |
| PERFORMANCE                | Seal strength average*           | cN/15mm | 330                                 | -                                     | 300                                  | -  |  |
|                            | Seal strength variation**        | %       | 4                                   | -                                     | 9                                    | -  |  |
|                            | Seal visual inspection           |         | Pass                                | Pass                                  | Pass                                 | Pass                                       |  |
|                            | Dye penetration (seal integrity) |         | Pass                                | Pass                                  | Pass                                 | Pass                                       |  |
|                            | Aseptic opening                  |         | Fiber free                          | Fiber free                            | Fiber free                           | Fiber free                                 |  |
|                            |                                  |         |                                     |                                       |                                      |  |  |

\* Sample size ; 60 blisters - 180° Peel angle. \*\*Average per pouch. \*\*\* Irradiation: 28.6 kGy. \*\*\* EO: Exposure Time 6h, Temperature 35°C at 50 % RH, EO concentration 1 kg/m<sup>3</sup>.

More extensive trials are currently being run at customer's facilities in Asia, Europe and United States of America.



#### NEXT STEPS

**POLYBOND AGP 90 gsm** is now currently under performance assessment and final testing review with several medical device manufacturers in Asia, Europe & United States of America.

The purpose of these testings is to gather more data on product behavior

- with different grades of PVC & APET films,
- different geometry of packaging
- different types of production lines (technology, age, maintenance policy)
- more significant statistical data

**POLYBOND BHP 85 GSM**, designed for sealing to HIPS material, is currently going through the first test assessment of performance as per POLYBOND AGP 90 gsm qualification protocol.

Interested in helping us learning more about product behavior, and qualifying it on your packaging lines?

#### STERIMED INFORMATION

From STERIMED, your Infection Control Partner:

- Largest range of sterilization packaging substrates in the world
- Sales in more than 70 countries
- Packaging education programmes for medical device industry and hospital staff

Accreditations:

• ISO 9001:2008 • ISO 14001:2004 • ISO 50001:2011 • ISO 13485:2003



Interested in helping us learning more about product behavior, and qualifying it on your packaging lines?

Send us an E-mail: contact@sterimed.us

ETHYPEL<sup>®</sup>, A BRAND OF



() ETHNEL () PROPNEL () STERISHEET

**STERIMED MARKETING DEPARTMENT** 32 avenue Pierre Grenier - 92100 Boulogne-Billancourt - FRANCE Phone: +33 157 75 92 75 - Fax: +33 1 57 75 92 93 - contact@sterimed.fr

www.sterimed.fr



STERIMED®, ETHYPEL®, PROPYPEL®, STERISHEET®, POLYBOND<sup>™</sup>, ETHYFLEX<sup>™</sup>, ETHYFORM<sup>™</sup>, CRYPTOPAK<sup>™</sup> STERIWEB®, 斯迪迈<sup>™</sup>, 舒卫安<sup>™</sup>, 蒸洁安<sup>™</sup>, Стеримед<sup>™</sup>, Этипель<sup>™</sup>, ПРОПИПЕЛЬ<sup>™</sup> are all brands of STERIMED

55

